# High-Q-Bond Bracket™

## LIGHT CURE BRACKET ADHESIVE SYSTEM

EN

### PROPERTIES

Bonds metal and ceramic brackets to tooth surfaces. Uses light cure adhesive technology to provide you with additional working time to ensure accurate bracket placement.

Small particle glass filled material designed specifically for orthodontic use that allows maximum penetration of filler into the finest screen mesh base. Releases Fluoride.

Contains the adhesive promoter 4-META, for strong reliable chemical bonding with tooth, metal and porcelain surfaces.

### INSTRUCTIONS FOR USE

1. Clean teeth with plain pumice using a rotating bristle brush or rubber cup.
2. Rinse and air-dry.
3. Etch for 20 seconds using **Q-Etch** or any other 36-38% phosphoric acid etching gel. Rinse and air-dry.
4. Isolate the bonded surfaces.
5. Apply a thin uniform layer of **High-Q-Bond Primer** on the etched enamel surface. **Gently spread the Primer by use of air!! Do Not rinse!** If desired, using a conventional halogen curing light or an LED light, polymerize the **Primer**. Curing the Primer is not mandatory. The primed surfaces should have a glossy appearance.
6. **NOTE:** Do not remove Primer with the high power suction!
7. Apply small amount of **High-Q-Bond Adhesive** on bracket base. **Adhesive** should be added to one bracket at a time.
8. Press bracket to tooth with a slight rotating motion. Final bracket adjustment must be made before **Adhesive** begins to cure from the ambient light (approximately 60 sec.).
9. Light cure entire bonded tooth surface for 20 seconds from two opposite directions.
10. When bonding to Porcelain, the bonded Porcelain surface should be etched for 60 seconds with Porcelain Etch, rinsed, dried thoroughly and primed with a Porcelain Silane. Please follow the instructions in the Porcelain Fix kit available separately. After removing residual Silane with a vigorous blast of oil free air, continue with step 5.

### CAUTION!

Refrigeration will increase viscosity, making syringing more difficult, and will retard setting time.

Contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible patients.

Wash thoroughly with soap and water after contact.

If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

### ESSENTIAL INGREDIENTS

- Urethan dimethacrylate.
- Triethylenglycoldimethacrylate.
- Multifunctional monomers.
- Photoinitiators.
- Fumed Silica.
- Adhesive accelerator.
- Ethanol (in Primer).

### STORAGE

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- Tightly recap after use
- To obtain SDS visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

**B.J.M. Laboratories Ltd.**  
12 Hassadna St., Industrial Park,  
Or-Yehuda 6022011, Israel  
tel. 972-74-7000111 fax. 972-3-7353020  
Info@bjmlabs.com [www.bjmlabs.com](http://www.bjmlabs.com)

**QUASER S.r.l.**  
Via Pian due Torri 19, 00146 Rome, Italy  
tel. 39-06-55282325 fax. 39-06-55290821  
info@quasersrl.it [www.quasersrl.it](http://www.quasersrl.it)

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