# Q-Crown™

## EN TEMPORARY CROWN & BRIDGE MATERIAL

Controlled radical self polymerizing, paste-paste composite, for fabrication of temporary crowns and bridges.

### INDICATIONS

Q-Crown is used to fabricate temporary crowns, bridges, inlays and onlays used in indirect restorative procedures.

### CONTRAINDICATIONS

Non known.

### PROPERTIES

- High compression resistance: 180-220 MPa (after 24h)
- High flexural strength 40-60 MPa (after 10 minutes), 120-150 MPa (after 24h)
- Great tensile and fracture strength
- Easy to remove during its elastic phase
- Does not heat up during polymerization
- Great color stability
- Easy to finish and polish

### INSTRUCTIONS FOR USE AUTO-MIX SYRINGE

- Bleed the syringe before installation of the automix

dispensing tip. • Mount a new, automix dispensing tip onto the Q-Crown dual-barrel syringe and bleed the syringe again. This insures complete flow of base and catalyst components. • Remove and discard the used dispensing tips and make sure to recap the syringe after each use. ☒

### PROPERTIES

Temporary crown and bridge fabrication can be made either directly or indirectly with impression trays. • The Q-Crown material reaches its elastic phase 60-90 seconds after start of mix, i.e. about 30-60 seconds after insertion in mouth. The impression must be removed from the mouth during this phase otherwise it may become difficult to remove. The working times at room temperature 23°C and 50% relative humidity, are as follows: • Initial setting in mouth (elastic phase) 1-1.5 minutes from start of mix. - During that period the material should be removed from the mouth. • Moderate setting occurs 2-2.5 minutes from start of mix. - After moderate setting remove the material from the impression tray, remove the thin inhibited layer

with alcohol (this makes the finishing easier) and carefully remove any excess or rough material in areas close to the gums. • Finishing and polishing are done only after the material has completely set (about 4-5 minutes from start of mix). • Cement the temporary using a non eugenol temporary cement (e.g. Q-Temp). • Brush a thin layer of the included Q-Crown Glaze on to the temporary and light cure each surface for 20 seconds.

### PRECAUTIONS

Q-Crown contains polymerizable monomers. Avoid prolonged or repeated contact with skin, oral soft tissues and eyes. Avoid prolonged inhalation. Do not take internally. Q-Crown should not be used with patients who have a history of severe allergic reaction to acrylate or methacrylate. Any compromise in oral hygiene or use of some mouth washes may lead to discoloration. Before use, lubricate the contacting hard surfaces to prevent breakage when the temporary is removed. Composite or Glass Ionomer cores must always be lubricated (e.g. with glycerine or Vaseline) to prevent them from adhering

to the Q-Crown material. Expiration date: 2 years when stored within the recommended temperatures.

### ESSENTIAL INGREDIENTS

#### Base & Catlyst:

- Triethylenglycol dimethacrylate.
- Photoinitiators.
- Multifunctional monomers.
- Polymerization activators.
- Fluoride releasing filler.
- Glass filler.
- Fumed Silica.

### STORAGE

- Do not store above 24°C (75°F).
- To obtain SDS visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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