

Q-Seal™

EN

LIGHT CURED, PIT & FISSURE SEALANT WITH FLUORIDE RELEASE


INDICATIONS

For the preventive sealing of pits and fissures.


PROPERTIES

- Easy & precise application by using direct syringe technique.
- Light-curing.
- Low viscosity, excellent flow properties, optimal wetting of fissures.
- Fluoride release.

INSTRUCTIONS FOR USE




1. Clean the enamel surface thoroughly with either cleansing paste containing neither oil nor fluoride, or with pumice and brush




5. Dry with air free of water and oil. The etched enamel should have a mat white appearance. Avoid contamination of the etched surface with saliva.



2. Isolate the working field, preferably with a rubber dam.




6. Mount and secure lock a new dispensing tip on the Q-Seal syringe and carefully syringe an adequate amount of Q-Seal into the fissures. Let the sealant flow until the entire fissure is covered. You can use a disposable brush, microbrush or the dispensing tip itself to spread the sealant and eliminate air bubbles.



3. Apply the etching gel, and let it react for 30 to 60 seconds.



4. Rinse thoroughly.



7. Cure sealant for 20 sec.

8. Check occlusion.

Q-Seal™

EN

CONTRAINDICATION

Sealing of pits and fissures with Q-Seal is contradicted:

- If pits or fissures are already affected with caries.
- If a dry working field cannot be established during the application and curing of Q-Seal.
- If the patient is known to be allergic to any of the ingredients of Q-Seal.

CAUTION!

Q-Seal contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact.

If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use

ESSENTIAL INGREDIENTS

- Urethan dimethacrylate.
- Multifunctional monomers.
- Triethylenglycoldimethacrylate.
- 4-Methacryloxyethyltrimellitic anhydride.
- Photoinitiators.
- Fumed Silica.

STORAGE

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

CE
0425



EC REP

B.J.M. Laboratories Ltd.
12 Hassadna St., Industrial Park,
Or-Yehuda 6022011, Israel
tel. 972-74-7000111 fax. 972-3-7353020
Info@bjmlabs.com www.bjmlabs.com

QUASER S.r.l.
Via Pian due Torri 19, 00146 Rome, Italy
tel. 39-06-55282325 fax. 39-06-55290821
info@quasersrl.it www.quasersrl.it

BJM LAB

2CRB005 Rev.5 10.2018