PRIMA 2000™

EN SINGLE STEP BONDING AGENT

INDICATIONS

- Direct composite restorations.
- Adhesive cementation of crown & bridges, including inlays and onlays.

PROPERTIES

• Prima 2000 is a dentin/enamel bonding agent which combines primer and adhesive into a single component for ease of application. • Prima 2000 is compatible with all common dental etching materials and forms an excellent bond with dentin and enamel, any self-cure or visible light-cured Composite, Compomer restorative material and with Composite Resin Cements. • Recent studies demonstrate that Prima 2000 is capable of delivering more than 20 MPa shear bond strength to dentin. • Scanning electron micrographs indicate that Prima 2000 application forms the essential hybrid zone on dentin.

INSTRUCTIONS FOR USE

Tooth Preparation:

Prepare the tooth according to clinical requirements.

Etching:

1. Etch enamel and/or dentin for 5-10 seconds. When bonding only to enamel, etch enamel for 20-30 seconds. ("Q-ETCH")
2. Rinse with a spray of oil-free water for 5-10 seconds.
3. Blot dry excess moisture from the surface. Verify that the surface has no visible moisture.

Application of Prima 2000:

1. Dispense Prima 2000 directly onto a brush or a microbrush by holding the brush tip directly over the nozzle of the tilted bottle.
2. Apply ample amount of Prima 2000, leave undisturbed for 15 sec. Remove solvents with a vigorous blast of air for at least 15 seconds or

until no visible movement of the resin is noted. If the surface appears incompletely coated, apply additional material.

3. Cure Prima 2000 for 20 seconds with a visible light-cure unit. (When the restoration retention area lies primarily within dentin, apply a second layer. Remove solvents immediately and light cure for 20 seconds).
4. Restore with visible light cured or self-cure composite resin of your choice or cement restoration.

CAUTION!

Prima 2000 requires a calcium hydroxide liner/base when working in close proximity to the pulp! Prima 2000 contains polymerizable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

ESSENTIAL INGREDIENTS

- Multifunctional methacrylates.
- Multifunctional monomers.
- Urethan dimethacrylate.
- Photoinitiators.
- Acetone.
- Ethanol.
- Adhesive accelerator.
- Colloidal Glass.

STORAGE

- Do not store above 24°C (75°F).
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- Tightly recap after use.
- To obtain SDS visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

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