# PORCELAIN FIX™

GB

## INDICATIONS:

BJM Porcelain Fix is to be used any time resin bonding to porcelain is required. It is used in:

- Preparation of porcelain crowns, veneers and inlays prior to cementation.
- Preparation of fractured porcelain or porcelain fused to metal crowns for bonding.

## PROPERTIES:

BJM Porcelain Fix is a complete kit for preparing porcelain veneers, inlays and fractured porcelain teeth for bonding. It consists of Porcelain Etch and Porcelain Silane. Porcelain Etch is a specially buffered, viscous, hydrofluoric acid that will etch porcelains of all type to produce a microscopic porous surface, which will provide a strong mechanical interlock with composite resin materials. Porcelain Silane is a single component Silane that enhances bonding values and durability when applied to etched porcelain.

## PREPARING PORCELAIN SURFACE FOR BONDING INSTRUCTIONS FOR USE:

1. In case of Intraoral procedure, isolate the area to be etched with a rubber dam. 2. When preparing fractured porcelain crowns for bonding, abrade the surface of the porcelain to be etched with a medium grit bur to remove the glazed surface and to make space for a sufficient layer of composite resin, for strength and aesthetics. Rinse all debris from the porcelain surface with copious amounts of water and dry with moisture and oil free compressed air. 3. Syringe the Porcelain Etch onto the porcelain surface. Leave undisturbed for 60 seconds. 4. Rinse the porcelain thoroughly with water to completely remove the etchant and dry well. The etched porcelain should have the frosted white appearance, if not, repeat the etching procedure. 5. Brush Porcelain Silane on the etched porcelain surface. Dry well with a stream of moisture and oil free compressed air.

Note: Porcelain bonding steps are addressed in the High-Q-Bond Adhesive Resin Cement illustrated in Instructions cards.

## SHELF LIFE:

Two years when stored below 24°C (75°F)

## CAUTION!

BJM Porcelain Etch must be used carefully and cautiously. Avoid contact with oral soft tissues, skin and eyes. If accidental contact occurs, flush with copious amounts of water and seek medical attention immediately.

In case of Intraoral procedure, we strongly recommend to Isolate the area to be etched with a rubber dam.

## ESSENTIAL INGREDIENTS:

Porcelain Etch: hydrofluoric acid: 8.0% • water: 90.5% • xantan gum: 1.5%

Porcelain Silane: ethyl alcohol: 97.0% • glycidoxypropyltrim ethoxysilane: 3.0%

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**B.J.M. Laboratories Ltd.**