

Dual-Cure Core Compatibility to DBA using Self-Cured and Self-Etching Activators

B. Zalsman¹, A. Valdman¹, K. Lizenboim¹, A. Suvorov¹, I. Suvorov¹, W.A. McHale², H. Dodiuk-Kenig³

¹BJM Laboratories Ltd., Or-Yehuda, Israel; ²Premier Dental Products company, Plymouth Meeting, PA, USA ; ³Shenkar College of Engineering and Design, Ramat-Gan, Israel

Introduction

Successful adhesion to dentin of dual-cured composite core build-up materials has been a challenge when used with many visible light-cured (VLC) dentin bonding agents (DBA). Manufacturers developed 'self-cured activators' (SCA) to enhance the bond compatibility of their dentin bonding agent to dual-cured composites.

Objective

To evaluate the bond compatibility (Shear Bond Strength SBS) of dual-cured composite core build-up materials to several commercially available DBA with their SCA and an experimental self-etching bond enhancer (EXP).

Methods

- Flat bonding sites were prepared on the buccal surfaces of 24 freshly extracted bovine teeth by grinding the teeth on a water-cooled abrasive wheel.
- 5 dual-cured composite core build-up materials (**Table 1**) were bonded to bovine dentin, according to the manufacturer instructions, utilizing 5 VLC DBA (**Tables 2,3**) with their SCA (**Table 4**):
 - with and without pre-treatment with:
 - self-etching primers (SEP) (**Table 5**)
 - an experimental self-etching bond enhancer (EXP), BJM Laboratories Ltd. (**Table 5**).
- Shear bond strength (SBS) was measured in accordance with ISO/TS11405:2003.
- The experimental results were statistically analyzed (N=10) by ANOVA (p<0.05).

Table 1. Dual-Cured Composite Core Build-Up Materials

| Component | Brand Name | Manufacturer |
|-----------|------------|-------------------------|
| CCAF | CompCoreAF | Premier Dental Products |
| LC | LuxaCore | DMG |
| CP | Core Paste | DenMat |
| TC | Ti-Core | EDS |
| QC | Q-Core | BJM Laboratories Ltd. |

Table 2. 5th generation Visible Light-Cured (VLC) Dentin Bonding Agents (DBA)

| Component | Brand Name | Manufacturer |
|-----------|--------------------|-------------------------|
| IB | IntegraBond | Premier Dental Products |
| OSP | Optibond Solo Plus | Kerr |
| P&BNT | Prime & Bond NT | Dentsply |
| PB | Prima 2000 | BJM Laboratories Ltd. |

Table 3. 6th generation Self Etching Visible Light-Cured (VLC) Dentin Bonding Agents (DBA)

| Component | Brand Name | Manufacturer |
|------------|--------------------|-----------------------|
| OSP | Optibond Solo Plus | Kerr |
| XIII (A&B) | Xeno III | Dentsply |
| CNXB | Contax Bond | DMG |
| PQB | Prima Quick Bond | BJM Laboratories Ltd. |

Table 4. Self-Cured Activators (SCA)

| Component | Brand Name | Manufacturer |
|-----------|------------------------------|-------------------------|
| ACA | Premier Auto Cure Activator | Premier Dental Products |
| OSPA | Optibond Solo Plus Activator | Kerr |
| NTA | Prime & Bond NT Activator | Dentsply |
| BACA | Auto-Cure Activator | BJM Laboratories Ltd. |

Table 5. Self-Etching Primers (SEP)

| Component | Brand Name | Manufacturer |
|-----------|---|-----------------------|
| OSPPR | Optibond Solo Plus Primer | Kerr |
| CNXP | Contax Primer | DMG |
| PQP | Prima Quick Prime | BJM Laboratories Ltd. |
| EXP | Experimental Self-Etching Bond Enhancer | BJM Laboratories Ltd. |

Results

The SBS (MPa) test results are presented in the following figures and tables.






Table 6. Shear Bond Strength [MPa] of Different Dual Cure Core Materials Bonded to Bovine Dentin Utilizing 5th Generation VLC Dentin Bonding Agents with and without their Self-Cured Activators and Experimental Bond Enhancer

| | CompCoreAF | CorePaste | LuxaCore | Ti-Core |
|----------------|------------|-----------|----------|----------|
| ETCH->IB+ACA | 6.8±4.0 | 5.1±3.2 | 5.7±2.4 | 7.2±2.8 |
| EXP->IB | 19.4±3.1 | 21.3±5.4 | 22.5±7.2 | 22.2±4.2 |
| ETCH->OSP+OSPA | 6.9±2.6 | 7.5±3.3 | 8.8±3.5 | 7.6±3.7 |
| EXP->OSP | 17.1±2.5 | 18.4±3.8 | 19.1±4.2 | 19.8±3.2 |

Table 7. Shear Bond Strength [MPa] of Different Dual Cure Core Materials Bonded to Bovine Dentin Utilizing 6th Generation VLC Dentin Bonding Agents with and without their Self-Etching Primers and Experimental Bond Enhancer

| | CompCoreAF | Q-Core |
|------------------|------------|----------|
| OSPPR → OSP | 16.4±2.0 | 18.6±2.7 |
| EXP → OSP | 17.1±2.5 | 23.1±2.6 |
| PQP → PQB | 15.8±1.8 | 16.7±2.4 |
| EXP → PQB | 21.9±2.1 | 22.4±3.0 |
| XIII (A&B) | 11.7±1.3 | 12.9±1.6 |
| EXP → XIII (A&B) | 16.5±2.1 | 14.4±2.2 |
| CNXP → CNXB | 1.0±0.2 | 1.3±0.2 |
| EXP → CNXB | 13.6±1.7 | 11.8±0.8 |

Discussion

BJM'S PQP was modified in order to enhance the compatibility of its bonding agent PQB to DC core materials. The modified formulation was later on referred to as the experimental self-etching bond enhancer (EXP).

The following modifications were incorporated into the PQP:

- Reduced % of solvents
- Added relatively lower molecular weight dimethacrylates
- Added difunctional methacrylate monomers and surface-modified, synthetic SiO₂-nano-spheres of 20 nm
- Added proprietary initiators

Conclusions

- The modified PQP, code name EXP, enhanced the SBS of PQB to DC core materials.
- All 5 core materials bonded to the EXP pre-treated, un-etched, bovine dentin without SCA exhibited significantly higher SBS compared to the same materials bonded with SCA to etched dentin.
- The EXP treated surfaces for all 6 generation DBA, exhibited higher SBS than the original DBA primers.