BJM Root Canal Sealer™

TWO-PASTE EPOXY-AMINE RESIN ROOT CANAL SEALER

EN

INDICATIONS

Permanent obturations of root canals combined with root canal points.

PROPERTIES

BJM Root Canal Sealer (BJM RCS) is a two-paste epoxy-amine resin root canal sealer in a dual auto-mixing syringe with subsequent characteristics:

- Outstanding sealing ability
- High radiopacity
- Non cytotoxic
- Long term stability
- Moderate flexibility that prevents cracking of fully cured material
- Low shrinkage
- Good flow properties
- Complies with ISO 6876 standard

Delivery device: 5ml auto mixing syringe with mixing element and Intra-oral tip. 

CONTRAINDICATIONS

Hypersensitivity to epoxy resins and amino groups. Attention: do not to use in case of patient who has allergic sensitivity to epoxies and amines.

In case of contact with eyes or skin:

wash immediately with plenty of water. Any excess of material must be removed with cotton applicator and then washed with water followed with suction.

For external use only! If swallowed, drink plenty of water. If any sickness symptoms seek medical assistance.

PRECAUTIONS

- Under certain storage conditions, BJM RCS Amine Paste may slightly separate. It has been shown that this process does not negatively affect the performance of the mixed product.
- Acid may slow epoxy polymerization (setting process). Please make sure to wash the canal with copious amounts of water and dry well prior to use of BJM RCS.
- You may observe some brown discoloration in the Amine Paste of BJM RCS. This is normal behavior of the material and has no effect on its mechanical or antibacterial properties. After mixing the discoloration will disappear and blend into the dentin shade color of the mixed material.

INSTRUCTIONS FOR USE

Prepare the root canal, clean and dry before use of BJM RCS.
Remove the cap by turning it 90° counter-clockwise and pulling it off. Bleed the syringe prior to use.
Attach the mixing tip to the syringe and then turn the tip 90° clockwise. Auto mixing syringe

ensures proper homogeneous mixing of two pastes and allows direct application of the sealer into the root canal.

After the application, the mixing tip must be removed by rotating it 90° counter clockwise and pulling it off. Then the auto mixing syringe must be closed with the original cap.

APPLICATION

As a standard technique, BJM RCS is used in combination with gutta-percha points. For filling techniques where most of the canal is obturated by endodontic point material, apply only a light coating of BJM RCS to the canal walls.

LATERAL CODENSATION TECHNIQUE

For the Master-Point-Technique, select a gutta-percha point (or alternatively a paper point or a reamer) in the size of the last instrument used during apical preparation. Wet the canal walls with BJM RCS through a pumping or simultaneously rotating movement in a counter clockwise direction of the point/reamer. Alternatively, apply BJM RCS onto the tip of a Lentulo spiral. Advance the Lentulo spiral slowly to the apex running at very low speed. Avoid the formation of air bubbles in the material and overfilling of the canal. Withdraw Lentulo very slowly, still running at low speed. Dip point into BJM RCS and insert them into the canal with a slow pumping motion.

3 methods of BJM Root Canal Sealer's application in warm gutta-percha obturations

1. Vertical Compaction by Schilder and Machtou Methodology

- Put BJM RCS on Master Cone (taper .04 or .06) and insert it 0.5-1 mm shorter than Working Length (WL) by a pumping and brushing movement.
- Heated pluggers (about 150°C) softened and cutting gutta-percha and afterwards Cold Hand plugger gradually compact it until the apical 5 mm plug is formed
- Put BJM RCS on the canal walls and fill with injectable gutta-percha.

2. Continuous Wave Condensation by Buchanan

- Put BJM RCS on Master Cone (taper .04 or 06) and insert it 0.5-1 mm shorter than Working Length (WL) by a pumping and brushing movement .
- One continuous movement of heated plugger (150°C-200°C) cutting gutta-percha Cone 5 mm shorter WL

- Put BJM RCS on the canal walls and fill with injectable gutta-percha.

3. Combination Technique by Solomonov

- Put BJM RCS on Master Cone (taper .02) and insert it to Working Length (WL) by a pumping and brushing movement.
- Deep Lateral Compaction from WL to the 5 mm from the WL with Ni-Ti spreader No. 25 and accessory gutta-percha points.
- Heating of gutta-percha points inside the canal up to 100°C.
- Compaction of the heated gutta-percha by hand pluggers Machtou.
- Repeat the procedure until the dense obturation of the canal.

REMOVAL OF BJM RCS

The conventional technique of removing the gutta-percha is applicable when BJM RCS is used with gutta-percha points.

WORKING TIME

BJM RCS will set within 48 hours at 37°C.

CLEANING OF INSTRUMENTS

All the instruments that were in contact with BJM RCS should be cleaned immediately after use with alcohol or acetone.

ESSENTIAL INGREDIENTS

Amine & Epoxy

- Zirconium Oxide.
- Fumed Silica.
- Poly (oxyethylene, oxypropylene) diamine.
- Bisphenol-A (epichlorohydrin) and Epoxy resin.
- Quaternum ammonium macromolecule.

STORAGE

- Do not store above 24°C (75°F).
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

 **B.J.M. Laboratories Ltd.**
12 Hassadna St., Industrial Park,
Or-Yehuda 6022011, Israel
tel. 972-74-7000111 fax. 972-3-7353020
Info@bjmlabs.com www.bjmlabs.com

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