PerfectTM

DESENSITIZER

EN

PERFECT DESENSITIZER IS USED

- To eliminate pain in exposed cervical areas not requiring restoration.
- To alleviate or prevent dentinal sensitivity after preparation of teeth to receive direct and Indirect restorations.

PROPERTIES

Perfect Desensitizer is used for the treatment of hypersensitive dentine. It achieves its effects by precipitation of plasma proteins, which reduces dentinal permeability and occludes the peripheral dentinal tubules. This inhibits the flow of fluid through the tubules which is the cause of sensitivity.

INSTRUCTIONS FOR USE

TREATMENT OF HYPERSENSITIVE DENTINE

The treatment of hypersensitive dentinal surfaces is only indicated in cases of severe pain in reactions to thermal, osmotic or tactile stimuli. In cases of mild pain, alternative treatments such as desensitizing toothpastes should be used.

1. Clean the dentine (e.g. with cleansing pumice), under local anaesthesia if necessary, and then rinse off with water.
2. Protect mucous membranes by using a rubber dam. Make sure that Perfect Desensitizer only comes into contact with the area to be treated.
3. Apply a small amount of Perfect Desensitizer to the dentinal surfaces using pellets or brushes and leave undisturbed for 30–60 seconds.
4. Then dry the surface carefully by applying a stream of compressed air until the fluid film has disappeared and the surface is no longer shiny.
5. Rinse thoroughly with water.

TREATMENT OF DENTINE PRIOR TO ADHESIVE BONDING OR ADHESIVE CEMENTATION PROCEDURES

1. Protect mucous membranes by using a rubber dam.
2. Cover cavities and/or areas near the pulp with suitable cavity linings before applying the desensitizer.
3. Remove moisture from the prepared dentine by applying a gentle stream of compressed air.
4. Apply Perfect Desensitizer with a pellet and leave for 30–60 seconds. Make sure only the smallest possible amount required is applied and that it only comes into contact with the area to be treated.
5. Dry the dentinal surface carefully by applying a stream of compressed air.
6. Rinse off thoroughly with water and apply suction.
7. Continue with Dentin Bonding or cement of your choice and complete restoration..

CAUTIONS

Perfect Desensitizer should only be used by a dentist

- Harmful, contains glutaraldehyde, 2-hydroxyethylmethacrylate.
- It is a local irritant and has toxic effects.
- Appropriate precautions should always be taken before using Perfect Desensitizer.
- Only dispense the smallest possible amount necessary for the treatment.
- Harmful by inhalation.
- Irritating to respiratory system.
- Irritating to skin.
- Risk of serious damage to eyes.
- May cause sensitization by inhalation.
- May cause sensitization by skin contact.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Wear suitable gloves.
- Before use, put on the protective goggles and cover the patient's eyes to protect against splashes of material.
- Do not breathe vapour/spray.
- Avoid contact with skin.

Perfect Desensitizer must not be used

- If the necessary precautions cannot be taken or the specified application technique cannot be used
- If the patient is allergic to any of the ingredients contained in Perfect Desensitizer.
- If the pain is due to pulpitis.
- If the pulp is exposed or the dentine adjacent to the pulp is perforated. Irritation of the pulp may be observed in the presence of deep cavities without suitable pulp protection. Those areas of cavities adjacent to the pulp should therefore be protected by a suitable cavity lining.

ESSENTIAL INGREDIENTS

- Multifunctional monomers.
- Glutaraldehyde.
- Water.

STORAGE

- Do not store above 24°C (75°F).
- Keep out of the reach of children.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- Tightly recap after use.
- To obtain SDS visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.

B.J.M. Laboratories Ltd.
12 Hassadna St., Industrial Park,
Or-Yehuda 6022011, Israel
tel. 972-74-7000111 fax. 972-3-7353020
Info@bjmlabs.com www.bjmlabs.com

QUASER S.r.l.
Via Pian due Torri 19, 00146 Rome, Italy
tel. 39-06-55282325 fax. 39-06-55290821
info@quasersrl.it www.quasersrl.it

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