High-Q-Bond Band™

High-Q-Bond Band™ Light-cure band cement

For professional dental use only

EN INDICATIONS

- Direct bonding of metal orthodontic bands to tooth surfaces.
- Occlusion adjustments to prevent deep bite.
- Occlusal interference with a particular orthodontic device.

PROPERTIES

- High-Q-Bond Band bonds chemically and mechanically to tooth surfaces.
- High-Q-Bond Band is a light-cure material - compatible with all curing units available on the market.
- High-Q-Bond Band depth of cure is >4 mm.
- Glass filled material designed specifically for orthodontic use.
- Excellent flow for filling the gap between the band and the tooth.
- Blue visible color. To see shade guide visit www.bjmlabs.com.
- High-Q-Bond Band has been on the market since 2007 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CAUTION!

- Refrigeration will increase viscosity, making syringing more difficult, and will retard setting time.
- High-Q-Bond Band contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible patients. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operating light. Tightly recap the syringe after use.
- High-Q-Bond Band contains materials which: sensitive to light / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.

COMPOSITION

Essential ingredients	% wt.	Essential ingredients	% wt.
Urethane dimethacrylate	<20%	Photoinitiators	<1%
Triethylenglycoldimethacrylate	<20%	Glass filler	<60%
2-hydroxyethylmethacrylate	<5%	Fumed silica	<5%
Polymerization activators	<1%	Titanium oxide	<5%

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FILLER CONTENT

The filler system comprises glass fillers and modified silica:

60-70 weight-% or 42-50 volume-%.

The particle size of the fillers ranges: <6 µm.

INSTRUCTIONS FOR USE

1. Clean teeth with plain pumice using a rotating bristle brush or rubber cup.
2. Rinse and air-dry.
3. Isolate the bonded surfaces.
4. Apply a small amount of High-Q-Bond Band Cement onto the tooth-bearing surfaces of the band.
5. Immediately seat the band onto the tooth with a slight rotating motion. A small amount of cement should exude around the entire band periphery. Final adjustment must be made before the cement begins to cure from the ambient light (approximately 60 seconds).
6. Light-cure entire bonded tooth surface for 20 seconds from two opposite directions.
7. High-Q-Bond Band Cement can be easily wiped away with a cotton roll before it starts setting or can be removed with a scaler after it has set.
8. Tightly re-cap syringe after use.

	min:sec
Light curing in a light-wave range of 440 – 470 nm Compatible with all curing units available on the market	00:20

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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