Cem-Implant™

Cem-Implant™

Non-eugenol, self-cure resin cement for implant retained crowns

Professional dental use only.

Patent Population: Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and there are no exclusions for specific patient groups.

EN INDICATIONS

- Non-eugenol, acrylic-urethane polymer based temporary cement for implant luting.
- Long-term cementation for permanent implant retained restorations and long-term provisionals.

PROPERTIES

- Self-cure cement.
- ISO 4049 Classification: Type 2, Class 1.
- Rigid final set that ensures marginal seal and firm retention.
- Available in 2 shades: gingival and white. To see shade guide visit www.bjmlabs.com.
- Cem-Implant has been on the market since 2004 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CONTRAINDICATIONS

- Not to be used on exposed pulp.
- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Do not use in excess.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

CAUTION!

- Though rare, some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use, remove from patient, and refer patient to a physician.
- Contains materials which: irritating to eyes / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.
- Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. You may also contact our EU Authorised Representative: Obelis S.A. (contact details below).

COMPOSITION

Essential ingredients – BASE	% wt.	Essential ingredients – CATALYST	% wt.
Urethan dimethacrylate	<60%	Urethan dimethacrylate	<60%
Triethyleneglycoldimethacrylate	<20%	Triethyleneglycoldimethacrylate	<40%
2-Hydroxyethylmethacrylate	<20%	Benzoyl Peroxide	<5%
Polymerization activators	<5%	Fumed Silica	<5%
Fumed silica	<5%	Zirconium oxide	<10%
Titanium oxide	<1%		
Zirconium oxide	<10%		

FILLER CONTENT

The filler system comprises modified silica and minerals:
10-20 weight-% or 2-10 volume-%.

The particle size of the fillers range: <6 µm.

INSTRUCTIONS FOR USE

- Only a small bead of Cem-Implant around the margin is necessary to securely hold your crown. This is especially true with highly retentive restorations. Using less cement saves money and allows for easy retrievals!
- Tips are for single patient use only to prevent cross-contamination between patients.
- Use one Detachable Manual Adapter per one syringe lifetime.

- Check fit, occlusion, and proximal contacts. Complete all finishing procedures.
- Dry and isolate the prepared abutments. Absolute dryness is not necessary, but obvious contamination should be avoided.
- Make sure the interior of the restoration is clean and dry. To facilitate future removal of implant restorations, lubricate the surface of the abutment with a thin layer of K-Y Jelly™ or other water-based lubricant.

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

4. Automix (When using a dual-barrel syringe with automix syringe mix tip): a) Bleed the syringe before installation of the automix syringe mix tip. Mount a new, automix syringe mix tip onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. b) Insert a small bead of Cem-Implant onto a microbrush or spatula.

Hand Mix (When using a dual-barrel syringe with Detachable Manual Adapter): a) Bleed the syringe before installation of the Manual Adapter. Mount a new manual adapter onto the dual-barrel syringe and bleed the syringe again. This ensures complete flow of Base and Catalyst components. b) Using a clean spatula or mixing stick, mix the two components into a homogenous paste.

Hand Mix (When using a two syringes system of Base and Catalyst): a) Extrude equal amount of Base and Catalyst onto the selected mixing surface. b) Using a clean spatula or mixing stick, mix the two components into a homogenous paste.

5. Apply a thin layer of cement along the inside wall of the cervical third of the restoration, making sure to completely coat the margin with cement. This approach helps to reduce the amount of excess cement expressed apically while providing the desired retrievability.

6. Firmly set the crown and either hold in place or have the patient gently close on a cotton roll until the initial set. The initial set is characterized by a distinct gel state. Excess cement is best removed in its gel state, after reaching the material-working time, with a scaler or explorer.

7. Tightly re-cap syringes after use.

8. For the material time indications, see table below.

Time indications:

The working and setting times depend on the ambient and oral temperature.

The times shown are based on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00

- Mixing options:
Note: Working and setting times complies to the mentioned above range without dependence on mixing options.

Mixing options	Ratio
Automix – 1 Automix syringe	1:1 – Equal amount of Base and Catalyst
Hand Mix – 1 Automix syringe with Detachable Manual Adapter	1:1 – Equal amount of Base and Catalyst
Hand mix - 2 separate syringes	1:1 – Equal amount of Base and Catalyst

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F)
- Keep product in original packaging
- For optimum freshness, keep refrigerated
- Do not freeze
- If refrigerated, allow the syringe to reach room temperature
- To obtain SDS or IFU visit www.bjmlabs.com
- Do not use after expiration date. See expiry date on the label
- Shelf life is 2 years after the date of manufacturing
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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