# High-Q-Bond Bracket™

## High-Q-Bond Bracket™

### Light-cure bracket adhesive system

*For professional dental use only.*

**Patient Population:** Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and there are no exclusions for specific patient groups.

#### EN INDICATIONS

- Direct bonding of metal and ceramic orthodontic brackets to tooth surfaces.

#### PROPERTIES

- High-Q-Bond Bracket bonds chemically and mechanically for a superior bond, while flowing into the tightest mesh without bracket drift.
- High-Q-Bond Bracket is a light-cure material - compatible with all curing units available on the market.
- High-Q-Bond Bracket depth of cure is >4 mm.
- Glass filled material designed specifically for orthodontic use that allows maximum penetration of filler into the finest screen mesh base.
- The viscosity of High-Q-Bond Bracket was designed to prevent adhesive run-on and bracket skating.
- High-Q-Bond Bracket has been on the market since 2007 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

#### CONTRAINDICATIONS

- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

#### CAUTION!

- Refrigeration will increase viscosity, making syringing more difficult and will retard setting time.
- High-Q-Bond Bracket contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible patients. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operating light. Tightly recap the syringe after use.
- High-Q-Bond Bracket Adhesive contains materials which: sensitive to light / irritating to eyes/ harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.
- Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/ or patient is established. You may also contact our EU Authorised Representative: Obelis S.A. (contact details below).

#### COMPOSITION – Essential ingredients

HIGH-Q-BOND BRACKET PRIMER	% wt.	HIGH-Q-BOND BRACKET ADHESIVE	% wt.	Q-ETCH ORTHO	% wt.
Urethane dimethacrylate	<70%	Urethane dimethacrylate	<20%	Phosphoric acid	37%
Triethyleneglycoldimethacrylate	<30%	Triethyleneglycoldimethacrylate	<10%	Water	<60%
Polymerization activatos	<1%	2-hydroxyethylmethacrylate	<5%	Thickening agent	<10%
Photoinitiators	<1%	Polymerization activators	<1%	Pigment	<0.1%
Ethanol	<20%	Photoinitiators	<1%		
		Glass filler	<70%		

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#### zHigh-Q-Bond Bracket Primer

The product does not contain filler.

#### High-Q-Bond Bracket Adhesive

The filler system comprises glass fillers and modified silica:

- 70-85 weight-% or 50-60 volume-%.
- The particle size of the fillers ranges: <6 um.

#### Q-Etch Ortho

The product does not contain filler.

#### INSTRUCTIONS FOR USE

##### Nota:

**If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time.**

- Clean teeth with plain pumice using a rotating bristle brush or rubber cup.
  - Rinse and air-dry.
  - Etch for 20 seconds using Q-Etch Ortho or any other 36-38% phosphoric acid etching gel. Rinse and air-dry.
  - Isolate bonded surfaces.
  - Apply a thin uniform layer of High-Q-Bond Primer on the etched enamel surface. Gently spread the Primer by use of air – Do Not Rinse! If desired, using a conventional halogen curing light or an LED light, polymerize the Primer. Curing the Primer is not mandatory. The primed surfaces should have a glossy appearance.
- Note:** Do not remove Primer with the high-power suction!
- Apply small amount of High-Q-Bond Adhesive on bracket base. Adhesive should be added to one bracket at a time.
  - Press bracket to tooth with a slight rotating motion. Final bracket adjustment must be made before Adhesive begins to cure from the ambient light (approximately 60 seconds).
  - Light-cure entire bonded tooth surface for 20 seconds from two opposite directions.
  - When bonding to Porcelain, the bonded Porcelain surface should be etched for 60 seconds with Porcelain Etch, rinsed, dried thoroughly, and primed with a Porcelain Silane. Please follow the instructions in the Porcelain Fix kit available separately. After removing residual Silane with a vigorous blast of oil free air, continue with step 5.
  - Tightly re-cap syringe after use.

	min:sec
Light curing in a light-wave range of 440 – 470 nm Compatible with all curing units available on the market	00:20

#### STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit [www.bjmlabs.com](http://www.bjmlabs.com).
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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