High-Q-Bond™

Dual-cure permanent adhesive resin cement

For professional dental use only.

Patient Population: Any human undergoing dental treatment related to the procedures listed below under "Indications". The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional and there are no exclusions for specific patient groups.

EN INDICATIONS

- Final permanent cementation of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns, and Maryland bridges.
- Metal alloys adhesive opaquer.
- Porcelain repair.
- Proven Color Stability over time.

PROPERTIES

- High-Q-Bond Adhesive Resin Cement is a multi-purpose, dual-cure, radiopaque, permanent adhesive resin cement.
- High-Q-Bond Adhesive Resin Cement classification according to ISO 4049 is Type 2, Class 3.
- High-Q-Bond has a radiopacity value above 2.5 mm Aluminium. (Radiopacity is defined as the ability of a material to be visible in x-ray photographs).
- High-Q-Bond SE forms ideal bond with metal alloys, porcelain, ceramic, dentin, enamel, amalgam.
- High-Q-Bond Adhesive Resin Cement has a chemical affinity for metal.
- High-Q-Bond Adhesive Resin Cement is a dual-cure cement – it has a self curing mechanism which can also be light cured for an immediate set and extra strength at the margins.
- Light curing High-Q-Bond Adhesive Resin Cement gives your restorations immediate stability and easy clean up of excess cement.
- High-Q-Bond features high retention and total margin integrity.
- High-Q-Bond Adhesive Resin Cement has been on the market since 2000 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

CONTRAINDICATIONS

- The use of the medical device is contraindicated if the patient is known to be allergic or sensitive to any of its ingredients.
- Not to be used on exposed pulp.
- Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of resin-based cements.

CAUTION!

- High-Q-Bond Adhesive Resin Cement contains polymerizable monomers which may cause allergic contact dermatitis in susceptible patients and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact. Wash thoroughly with soap and water after contact. If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.
- High-Q-Bond Adhesive Resin Cement contains materials which: sensitive to light / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact. For complete symbols guide see below.
- Any serious incident that occurs in relation to this device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/ or patient is established. You may also contact our EU Authorised Representative: Obelis S.A. (contact details below).

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Essential ingredients BASE	% wt.	Essential ingredients CATALYST	% wt.
Polymerization activators	<1%	Triethylenglycoldimethacrylate	<20%
Photoinitiators	<1%	4-Methacryloxyethyltrimellitic anhydride	<5%
Urethan dimethacrylate	<40%	Benzoyl Peroxide	<1%
Glass filler	<50%	Glass filler	<60%
Fumed silica	<5%	Fumed Silica	<5%
Titanium oxide	<1%	Titanium oxide	<1%

FILLER CONTENT

The filler system comprises glass fillers and modified silica:

50-65 weight-% or 35-47 volume-%.

The particle size of the fillers range: <6 um.

INSTRUCTIONS FOR USE

Nota:

- Use are for single patient use only to prevent cross-contamination between patients.**
- Tips one Detachable Manual Adapter per one syringe lifetime.**

- Depending on the clinical requirement use the corresponding illustrated card.
- Time indications:
 - The working and setting times depend on the ambient and oral temperature.
 - The times shown are based on conditions relevant for practice.

	min:sec
Working time at ambient temperature	>01:30
Intra-oral setting time	<05:00
Light curing in a light-wave range of 440 – 470 nm <ul style="list-style-type: none">Compatible with all curing units available on the market	00:20

Mixing options:

Nota: working and setting times complies to the mentioned above range without dependence on mixing options.

Mixing options	Ratio
Automix – 1 Automix syringe	1:1 – Equal amount of Base and Catalyst
Hand Mix – 1 Automix syringe with Detachable Manual Adapter	1:1 – Equal amount of Base and Catalyst

STORAGE AND DISPOSAL

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- If refrigerated, allow the syringe to reach room temperature.
- To obtain SDS or IFU visit www.bjmlabs.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

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