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Dr. Michael Solomonov
DMD, Endodontist

**Department of Endodontics
Sheba Hospital, Tel-Hashomer, IDF Israel**

Dr. Michael Solomonov graduated from Tel Aviv University School of Dental Medicine in 1994. Dr. Solomonov received post-graduate endodontic training at the Endodontic department of Hebrew University, Hadassah School of Dental Medicine, in Jerusalem, Israel from 1999 to 2002 and served there as clinical instructor till 2010. He is Israel board certified endodontist from 2003 and serves as Examiner at Israel Board of Endodontics of the Scientific Council of the Israel Dental Association from 2008.

He retains a private endodontic specialist practice in Tel Aviv.

Dr. Solomonov is actively lecturing on endodontics worldwide. Dr. Solomonov is currently director of postgraduate endodontic program in Department of Endodontics, Sheba Hospital, Tel-Hashomer, Israel.

In my opinion Insoluble macromolecules as prevention of biofilm can lead us to a new level to solve the problem of coronal seal.

Dr. Michael Solomonov
DMD, Endodontist, Israel

BJM LAB

B.J.M. Laboratories Ltd.
12 Hassadna St., Industrial Park, Or-Yehuda 6022011,
ISRAEL Fax. 972-3-7353020 www.bjmlabs.com

Diadent, 13 Rue Jean Lolive, Z.I. La Fosse
A La Barbriere, AULNAY 93600, FRANCE
Tel. 331-48012710, Fax. 331-42465230

2ABR031BJM Rev.1 02.2017

BJM ROOT CANAL SEALER™

Two-paste epoxy-amine resin root canal sealer

INDICATIONS

- Obturation of root canals together with gutta-percha points.

PROPERTIES

- Extremely High radiopacity.
- Excellent wettability and flow properties.
- Outstanding sealing ability.
- IABT Antibacterial technology.
- Non-cytotoxic.
- Long-term stability.
- Moderate flexibility that prevents cracking of fully cured material.
- Low shrinkage.
- Automix Syringe - Saves application time; guarantees consistent mix.

TECHNICAL DATA

BJM RCS will set within 48 hours at 37°C
Shelf Life 2 Years

SCIENTIFIC PAPERS

1. Antibacterial mechanism of novel endodontic sealer, D. Kesler Shvero, N. Zaltsman, E. Weiss, N. Beyth, Hadassah School of Dental Medicine, Hebrew University, IADR Israeli Division Meeting, Tel-Aviv, June 2013.
2. Root canal sealers as Biofilm prevention: facts and speculations, M. Solomonov, Эндодонтия, Том VII, No. 1-2, 2014.
3. Antibiofilm Activity of Epoxy Sealer with Quaternary Ammonium Macromolecule, T. Becker, M. Solomonov, N. Sterer, R. Bar-Ness, A. Levin, A. Shemesh, The Maurice and Gabriela Goldschleger School of Dental Medicine Tel Aviv University, Program Number 0219, PER-IADR Congress, Jerusalem, Israel, 2016.

PACKAGING & ORDER INFORMATION


- **Item # 400200**
- 1 Automix Syringe 5 ml of BJM RCS
- 10 Automix Syringe Mix Tips and Inta Oral Tips
- 1 Mixing Pad

IABT incorporation into dental polymers prevents bacterial growth and biofilm formation.


BJM LAB

Case 1.

Complex anatomy and healing


Chronic Apical Abscess, Two visit treatment finished after Sinus Tract disappearance. 2 Mesial and 3 Distal canals prepared by Combination of Rotary NiTi and SAF' obturated by combination technique with guthapercha and BJM sealer. 6 month Follow up revealed healing lesion

Case 4.

Dense in Dente

Treatment of Dense in Dente type 3b with Big PERIAPICAL lesion after analyses of CBCT. Main oval-C-shape canal prepared by SAF system, Round Dense canal prepared by Rotary NiTi. Combination technique of obturation with prominent warm component in main canal by Guthapercha and BJM RCS

Case 2.

Internal Resorption

One visit treatment of internal resorption with SAF system. Combination technique of obturation: Lateral compaction with component of warm compaction: guthapercha and BJM sealer

Case 5.

Broken instrument removal

Bypass of broken instrument in oval ML canal with removal during canal instrumentation. Combination technique of obturation by guthapercha and BJM RCS

Case 3.

Rare Anatomy

3 separate Mesial canals and 2 distal instrumented by rotary NiTi ,Obturated by Lateral Compaction with gutta-percha and BJM RCS