

Bisphenol-A Free Alternatives For Dental Polymer-Based Materials


Kira Lizenboim¹, Hanna Dodiuk², Ira Suvorov¹, Barry Zalsman¹

¹*B.J.M. Laboratories Ltd., Or Yehuda, ISRAEL*

²*Shenkar College, The Pernick Faculty of Engineering, Ramat Gan, ISRAEL*

Rationale:

- **Bisphenol-A (BPA)** is suspected to be an endocrine disrupter (resembling Estrogen hormone)
- Current polymeric dental materials are based on **BPA derivatives**, e.g. **Bisphenol-A Diglycidylether Methacrylate (Bis-GMA)** which may leach out unreacted monomers and its degradation products
- The growing international concern regarding the presence of BPA in commercial products has led to many studies of its effect on human health
- However, this topic is still controversial because regulators think that the methods used in many of those studies haven't been fully validated




Objectives:



- To evaluate **Bisphenol-A free** alternatives for potential use in dental polymer-based materials:
 - Use of synthetic commercial oligomers
 - Newly synthesized monomers from renewable resources

Introduction:

- Dental restorations
 - Repair a damaged tooth



- Key properties:
 - Biocompatibility
 - Excellent bonding to tooth structure
 - Endurance within the aggressive oral environment
 - Dimensional stability
 - Easy handling and Dentist friendly
 - Matching natural tooth color



- The dental restorative system
 - Resin (BPA free)
 - Initiators
 - Fillers
 - Coupling agents (Resin - Filler)



Methods:

Synthetic Alternatives:


- 13 commercial oligomers produced by different manufacturers were studied as resins for dental composites:
 - Ethoxylated (4) Bisphenol A Dimethacrylate
 - Difunctional Polyester Acrylate
 - **Difunctional Aliphatic Urethane Acrylates (DAUA)**
 - Difunctional Aliphatic Urethane Methacrylate
 - Difunctional Aliphatic Polyester Tri-Urethane Acrylate
 - Urethane Di-Methacrylate
 - Trifunctional Aliphatic Polyester Urethane Acrylate
 - Isocyanurate
 - Glycerol Dimethacrylate / Maleate
 - Urethane Dimethacrylate Elastomer
 - **Hexafunctional Aliphatic Urethane Acrylate (HAIUA)**
 - **Hexafunctional Aromatic Urethane Acrylate (HArUA)**

Renewable Resource Monomers:

- Bile Acid (Cholic Acid)
 - Synthesis of a tri-methacrylate- monomer of cholic acid was carried out and verified by different methods as **methacrylated derivative of cholic acid methyl ester (CAME)**
 - The cholic acid derivative (CAME) were evaluated as Bis-GMA replacement



- Commercially available vinyl macromer (**EcoMer**)



Physical, mechanical and adhesive properties of HArUA, HAIUA and one of the two tested DAUA excelled in all formulations.

Testing Methods:

Nuclear Magnetic Resonance (NMR), FTIR, DSC - verification of the synthesized cholic acid derivative;
Compressive, Flexural, Shear Bond Strength; Water Sorption, Light Curing Time - Degree of Conversion, Shrinkage, Flow Properties – verification of composites and sealants.
All data was statistically analyzed by the analysis of variance (ANOVA) method, at a significance level set at $p < 0.005$.

Preparation of Light Cure Composite and Sealant:

- 25 different **composite formulations** based on the five different oligomers:
 - Different ratios of Resin/Bis-GMA / TEGDMA (11-21)
 - Photoinitiator camphorquinone (CQ) and accelerator ethyl-4-dimethylaminobenzoate (EDB) were added to the composite mixtures
 - Addition of the photoinitiators TPO and Irgacure 819 to the composite mixtures was evaluated
 - The Filler comprised fumed silica (Cab-O-Sil) and barium-aluminum-borosilicate glass powder (Dental Glass) of 4 different grinding sizes (0.4, 1.5, 6 μm and 180 nm) (F1-F5)
 - Hydroxyapatite (HAp) addition (nanostructured microparticles having specific surface area $> 100 \text{ m}^2/\text{g}$) to the Filler mixture was evaluated (F6)
- The Filler content of the composites comprised of 3% untreated fumed silica (single particle size – 17nm, aggregate size 200-300nm) (Cab-O-Sil) and of 97% Dental Glass having a refractive index similar to the composite resin matrix (1.53) and treated with γ -Methacryloxypropyl-tri-methoxy-Silane. In some formulations 5% Dental Glass were replaced by 5% hydroxyapatite (HAp) (fine powder). Each formulation contained 75wt% of filler.
- Light cure dental composite restorative material ProFil, (made by Silmet Ltd.) was used as reference. ProFil contains Bis-GMA and TEGDMA resins and 75wt% glass filler.
- 3 different light cure **pit & fissure sealant formulations** based on three different oligomers containing photoinitiator camphorquinone (CQ) and accelerator ethyl-4-dimethylaminobenzoate (EDB).
- Each sealant formulation contained 3wt% of untreated fumed silica (single particle size – 17nm, aggregate size 200-300nm) (Cab-O-Sil) or 5% hydroxyapatite (HAp) (fine powder) (nanostructured microparticles having specific surface area $> 100 \text{ m}^2/\text{g}$) as a filler.
- Light cure pit & fissure sealant Q-Seal, (made by B.J.M. Laboratories Ltd.) was used as reference. Q-Seal contains Bis-GMA and TEGDMA resins and 3wt% filler.

Results:

Fillers Compositions

Filler Mixture Number	Micron Particle Sizes (wt%)			Nano Particle Sizes (wt%)		HAp (wt%)	Total	$\mu\text{m}/\text{nm}$ ratio
	Dental Glass (6 μm)	Dental Glass (1.5 μm)	Dental Glass (0.4 μm)	Dental Glass (180nm)	Cab-O-Sil (17 nm)			
F1	97	0	0	0	3	0	100	32
F2	55	28	0	14	3	0	100	5
F3	55	28	0	14	3	0	100	5
F4	55	28	0	14	3	0	100	5
F5	43.4	32.6	21.6	0	2.4	0	100	41
F6	41.0	30.0	21.6	0	2.4	5.0	100	41

Mechanical Properties of Selected Oligomers / Fillers Light Cure Composite Formulations

containing photoinitiator CQ and accelerator EDB

Mixture number	Oligomer mixture number	Oligomer mixture description	Filler mixture number	Compressive Strength, MPa	Flexural Strength, MPa
11-F1	11	HArUA:TEGDMA:Bis-GMA	F1	202±18 ^a	189±33 ^b
11-F2	11	HArUA:TEGDMA:Bis-GMA	F2	257±22 ^a	162±23 ^b
11-F3	11	HArUA:TEGDMA:Bis-GMA	F3	236±10 ^a	220±21 ^b
11-F4	11	HArUA:TEGDMA:Bis-GMA	F4	135±40 ^a	212±35 ^b
11-F5	11	HArUA:TEGDMA:Bis-GMA	F5	260±21 ^a	239±31 ^b
12-F1	12	EBAD:TEGDMA:Bis-GMA	F1	180±38 ^a	182±8 ^b
13-F1	13	HArUA:TEGDMA=70:30	F1	149±29 ^a	177±20 ^b
13-F5	13	HArUA:TEGDMA=70:30	F5	272±28 ^a	239±31 ^b
14-F1	14	HArUA:TEGDMA=80:20	F1	167±40 ^a	221±14 ^b
14-F4	14	HArUA:TEGDMA=80:20	F4	181±19 ^a	161±28 ^b
14-F5	14	HArUA:TEGDMA=80:20	F5	294±21 ^a	227±13 ^b
14-F6	14	HArUA:TEGDMA=80:20	F6	245±13 ^a	213±30 ^b
15-F1	15	CAME:TEGDMA=60:40	F1	70±8 ^a	110±11 ^b
15-F5	15	CAME:TEGDMA=60:40	F5	130±9 ^a	64±12 ^b
16-F5	16	DAUA:TEGDMA=80:20	F5	246±18 ^a	178±15 ^b
16-F6	16	DAUA:TEGDMA=80:20	F6	229±19 ^a	194±18 ^b
17-F5	17	HAIUA:TEGDMA=80:20	F5	220±18 ^a	217±22 ^b
17-F6	17	HAIUA:TEGDMA=80:20	F6	223±11 ^a	243±26 ^b
ProFil (Ref.)				223±33 ^a	274±52 ^b

containing photoinitiator TPO in addition to CQ and EDB

14-F5	14	HArUA:TEGDMA=80:20	F5	238±23 ^a	201±35 ^b
16-F5	16	DAUA:TEGDMA=80:20	F5	238±32 ^a	191±14 ^b
17-F5	17	HAIUA:TEGDMA=80:20	F5	210±22 ^a	210±22 ^b

Physical Properties of Investigated Bisphenol-A Free Dental Composites

Property	Standard Requirements	CAME based	HArUA based	HAIUA based	DAUA based
Biocompatibility	Material should show no cytotoxicity	N/E	N/E	N/E	N/E
Light curing time @ 23°C	Not more than 30 sec	20 sec	20 sec	20 sec	20 sec
Degree of Conversion, %	Degree of conversion should reach at least 50% 5 min after the start of irradiation by a visible light-curing lamp	44	45	76	40
Water sorption, $\mu\text{g}/\text{mm}^3$	Maximum 40 $\mu\text{g}/\text{mm}^3$	40	23	15	15
Solubility, $\mu\text{g}/\text{mm}^3$	Maximum 7.5 $\mu\text{g}/\text{mm}^3$	14	2	3	4
Flexural strength, MPa	Minimum 50 MPa	64	227	217	178
Compressive strength, MPa	Minimum 50 MPa	130	294	220	246
Radiopacity, % Al	Minimum 100 % Al	✓	✓	✓	✓
Adhesion (SBS to etched bovine enamel, MPa)	Minimum 15 MPa to etched enamel using a bonding agent	33	38	36	35
Volumetric shrinkage, %	Maximum 7 %	<3	<3	<3	<3
Wear resistance	Maximum volume loss of 10 μm per 39,000 3-body wear cycles	N/E	N/E	N/E	N/E
Handling property: viscosity vs. shear rate	While the viscosity decreases under an increasing shear rate, once the force is removed, the viscosity should return to nearly the higher viscosity	✓	✓	✓	✓
Shade (color) stability	Material should be visually homogeneous and match the defined color	✓	✓	✓	✓

Standard deviations are \pm (a) 10%, (b) 15% and (c) 20%. Values with the different superscript letters (a, b, c) in the same column are statistically different ($p < 0.005$).

N/E – Not Established in our laboratory.

containing photoinitiator Irgacure 819 in addition to CQ and EDB

14-F5	14	HArUA:TEGDMA=80:20	F5	227±38 ^a	199±13 ^b
16-F5	16	DAUA:TEGDMA=80:20	F5	192±21 ^a	171±18 ^b
17-F5	17	HAIUA:TEGDMA=80:20	F5	250±21 ^a	230±14 ^b

Adhesive and Mechanical Properties of Selected Oligomers / Filler Light Cure Pit & Fissure Sealant Formulations

Oligomer / Filler mixture description	Compressive Strength, MPa	Flexural Strength, MPa	Shear Bond Strength to Etched Bovine Enamel, Mpa / Failure Mode
HArUA:TEGDMA=80:20 / Cab-O-Sil	190±36 ^a	200±26 ^b	50±6 ^b / Cohesive
HArUA:TEGDMA=80:20 / HAp	228±19 ^a	185±14 ^a	54±9 ^b / Cohesive
DAUA:TEGDMA=80:20 / Cab-O-Sil	154±27 ^a	191±37 ^b	42±5 ^b / Mixed
DAUA:TEGDMA=80:20 / HAp	330±65 ^a	204±42 ^a	49±6 ^b / Cohesive
HAIUA:TEGDMA=80:20 / Cab-O-Sil	201±28 ^b	228±11 ^a	48±6 ^b / Cohesive
HAIUA:TEGDMA=80:20 / HAp	244±36 ^b	247±20 ^b	49±6 ^b / Cohesive
Q-Seal (Ref.)	191±54 ^a	205±33 ^a	45±5 ^b / Cohesive

Standard deviations are \pm (a) 10%, (b) 15% and (c) 20%. Values with the different superscript letters (a, b, c) in the same column are statistically different ($p < 0.005$).

Conclusions:

- The objective was to replace the Bis-GMA resin
- Maintain its key properties
- Changing as little as possible the composition (Resin/TEGDMA ratios, initiators and filler)
- The experimental results showed that Bis-GMA can be substituted by
 - Tri-methacrylated derivative of **cholic acid** methyl ester (referred to as **CAME**)
 - Aromatic-urethane hexafunctional oligomer (**HArUA**)
 - Hexafunctional Aliphatic Urethane Acrylate (**HAIUA**)
 - Difunctional Aliphatic Urethane Acrylates (**DAUA**)