BJM Root Canal Sealer™

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Michael Solomonov
DMD

Department of Endodontics
Sheba Hospital, Tel-Hashomer, Israel

Dr. Michael Solomonov graduated from Tel Aviv University School of Dental Medicine in 1994.

Dr. Solomonov received post-graduate endodontic training at the Endodontic department of Hebrew University, Hadassah School of Dental Medicine, in Jerusalem, Israel from 1999 to 2002 and served there as clinical instructor till 2010. He is Israel board certified endodontist from 2003 and serves as Examiner at Israel Board of Endodontics of the Scientific Council of the Israel Dental Association from 2008.

He retains a private endodontic specialist practice in Tel Aviv.

Dr. Solomonov is actively lecturing on endodontics worldwide. Dr. Solomonov is currently director of postgraduate endodontic program in Department of Endodontics, Sheba Hospital, Tel-Hashomer, Israel.

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- Hume W. R. Influence of dentine on the pulpward release of eugenol or acids from restorative materials. *Journal of Oral Rehabilitation* 21, no. 4 (1994): 469–473.
- Nair P. N. R., Stephane Henry, Victor Cano and Jorge Vera. Microbial status of apical root canal system of human mandibular first molars with primary apical periodontitis after «one-visit» endodontic treatment. *Oral Surgery, Oral Medicine, Oral Pathology, Oral Radiology, and Endodontology* 99, no. 2 (2005): 231–252.
- Saleh I. M., I. E. Ruyter, M. Naapasalo, and D. Orstavik. Survival of *Enterococcus faecalis* in infected dentinal tubules after root canal filling with different root canal sealers in vitro. *International Endodontic Journal* 37, no. 3 (2004): 193–198.
- Heling Ilana and Nicholas Paul Chandler. The antimicrobial effect within dentinal tubules of four root canal sealers. *Journal of endodontics* 22, no. 5 (1996): 257–259.
- Ray H. A., M. Trope. Periapical status of endodontically treated teeth in relation to the technical quality of the root filling and the coronal restoration. *International Endodontic Journal* 28, no. 1 (1995): 12–18.
- Magura Mark E., Abdel H. Kafrawy, Cecil E. Brown, Carl W. Newton. Human saliva coronal microleakage in obturated root canals: an in vitro study. *Journal of Endodontics* 17, no. 7 (1991): 324–331.
- Ricucci Domenico and Jose F. Siqueira Jr. Biofilms and apical periodontitis: study of prevalence and association with clinical and histopathologic findings. *Journal of Endodontics* 36, no. 8 (2010): 1277–1288.
- Kenawy E.R., Worley S.D., Broughton R. The chemistry and applications of antimicrobial polymers: a state-of-the art review. *Biomacromolecules*, 2007; 8(5):1359–1384.
- Kishen A., Shi Z., Shrestha A. et al. An inves-

tigation on the antibacterial and antibiofilm efficacy of cationic nanoparticles for root canal disinfection. *J Endod* 2008; 34:1515–20.

10. Upadya Megha, Annie Shrestha and Anil Kishen. Role of efflux pump inhibitors on the antibiofilm efficacy of calcium hydroxide, chitosan nanoparticles, and light-activated disinfection. *Journal of endodontics* 37, no. 10 (2011): 1422–1426.

11. Wu D., Fan W., Kishen A., Gutmann J.L., Fan B. Evaluation of the Antibacterial Efficacy of Silver Nanoparticles against *Enterococcus faecalis* Biofilm. *Journal of endodontics*, 2014; 40(2), 285–290.

12. Beyth N., Yudovin-Farber I., Bahir R., Domb A.J., Weiss E.I. Antibacterial activity of dental composites containing quaternary ammonium polyethyleneimine nanoparticles against *Streptococcus mutans*. *Biomaterials* 2006;27:3995–4002.

13. Yudovin-Farber Ira, Nurit Beyth Ervin I. Weiss and Abraham J. Domb. Antibacterial effect of composite resins containing quaternary ammonium polyethyleneimine nanoparticles. *Journal of Nanoparticle Research* 12, no. 2 (2010): 591–603.

14. Nisimov N. Zaltsman, D. Kesler, Ei. Weiss, N. Beyth. The antibacterial effect of a core resin buildup incorporating QPEII. Manuscript in preparation.

15. E. Varon-Shahar, N. Beyth. Antibacterial activity of the orthodontic cements incorporating polyethyleneimine against *Streptococcus Mutans*. Manuscript in preparation

16. Hashimoto Masanori, Hiroki Ohno, Hidehiko Sano, Franklin R. Tay, Masayuki Kaga, Yoshiyuki Kudou, Haruhisa Oguchi, Yoshima Araki and Minoru Kubota. Micromorphological changes in resin?dentin bonds after 1 year of water storage. *Journal of Biomedical Materials Research* 63, no. 3 (2002): 306–311.

17. Iris Slutzky-Goldberg, Hagay Slutzky, Michael Solomonov, Joshua Moshonov, Ervin I. Weiss and Shlomo Matalon. Antibacterial Properties of Four Endodontic Sealers. *Journal of Endodontics*; 34 (2008): 735–738.

18. Abramovitz Itzhak, Nurit Beyth, Yafit Paz, Ervin I. Weiss, and Shlomo Matalon. Antibacterial temporary restorative materials incorporating polyethyleneimine nanoparticles. *Quintessence international* 44, no. 3 (2012): 209–216.

19. D. Kesler Shvero, N. Zaltsman, E. Weiss, N. Beyth Antibacterial mechanism of novel endodontic sealer. Manuscript in preparation.

20. Barros J., Silva M.G., Rocas I.N., Goncalves L.S., Alves F.F., Lopes M.A. & Siqueira Jr, J.F. Antibiofilm effects of endodontic sealers containing quaternary ammonium polyethyleneimine nanoparticles. *Journal of Endodontics*, 2014; Aug; 40(8):1167–71.

21. Lockman P.R. et al. Nanoparticle surface charges alter blood-brain barrier integrity and permeability. *J Drug Target*. 2004; 12(9–10): 635–641.

22. FDA, 2010. Center for Drug Evaluation and Research MAPP. Office of Pharmaceutical Science. Reporting Format for Nanotechnology.

23. Rocks S.S., Pollard R.D., Levy L., Harrison P., Handy R. Comparison of risk assessment approaches for manufactured nanomaterials. 2008; Defra, London.

24. Chaudhry Q., Bouwmeester H. and Hertel R.F. (2010) The Current Risk Assessment Paradigm in Relation to the Regulation of Nanotechnologies, In: G.A. Hodge, D.M. Bowman and A.D. Maynard (eds), *International Handbook on Regulating Nanotechnologies*. Cheltenham: Edward Elgar, 124–143.

25. D'Antonio N.N., Rihs J.D., Stout J.E., Yu V.L. Computer keyboard covers impregnated with a novel antimicrobial polymer significantly reduce microbial contamination. *Am J Infect Control*. 2013 Apr; 41(4):337–9.

26. Shemesh A., A. Levin, Ben Itzhak, B. Katzenell, M. Solomonov. Comparisons of the physical properties of 3 epoxy resin-based root canal sealers – a novel one and two old. Manuscript in preparation.

Root Canal Sealers as Biofilm Prevention: Facts and Speculations

M. Solomonov

Abstract

The present paper is dedicated to the influence of endodontic sealer composition on the microbial biofilm within the root canal. The data of the recent studies concerning the use of insoluble macromolecules and nanoparticles as antibacterial agents included in the root canal sealers are also cited. Furthermore this article summarizes the information about its mechanisms of action on the endodontic infection, advantages of these new disinfection agents and the potential risks of its use.

Key words: root canal sealer, obturation of the root canal, biofilm, disinfection, nanoparticles.

BJM Root Canal Sealer™

TWO-PASTE EPOXY-AMINE
RESIN ROOT CANAL SEALER

Indications

- Obturation of root canals together with gutta-percha points.

Properties

- Extremely high radiopacity.
- Excellent wettability and flow properties.
- Outstanding sealing ability.
- IABT Antibacterial technology.
- Non-cytotoxic.
- Long-term stability.
- Moderate flexibility that prevents cracking of fully cured material.
- Low shrinkage.
- Automix Syringe - Saves application time; guarantees consistent mix.

Scientific Papers

1. Antibacterial mechanism of novel endodontic sealer, D. Kesler Shvero, N. Zaltsman, E. Weiss, N. Beyth, Hadassah School of Dental Medicine, Hebrew University, IADR Israeli Division Meeting, Tel-Aviv, June 2013.
2. Root canal sealers as Biofilm prevention: facts and speculations, M. Solomonov, **Эндодонтия**, Том VII, No. 1-2, 2014.
3. Antibiofilm Activity of Epoxy Sealer with Quaternary Ammonium Macromolecule, T. Becker,

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IABT incorporation into dental polymers prevents bacterial growth and biofilm formation.

M. Solomonov, N. Sterer, R. Bar-Ness, A. Levin, A. Shemesh, The Maurice and Gabriela Goldschleger School of Dental Medicine Tel Aviv University, Program Number 0219, PER-IADR Congress, Jerusalem, Israel, 2016.

4. Dr. Michael Solomonov, Clinical Cases Report, March 2017.
5. Evaluating the physical properties of one novel and two well-established epoxy resin-based root canal sealers, M. Solomonov, J.B. Itzhak, Quintessence Publishing Endo 2017; 11 (4): 285-290.
6. Antibiofilm activity of epoxy sealer incorporated with quaternary ammonium macromolecule, T. Becker, N. Sterer, R. Bar-Ness, T. Toledano and M. Solomonov Evidence-Based Endodontics 2019; 4 (1): 1-6.
7. Dr. Gabriel Kaplan, Clinical Cases Report, March 2019.

Technical Data

BJM RCS will set within 48 hours at 37°C
Shelf Life 2 Years

Packaging & Order Information

- Item # 400200
- 1 Automix Syringe 5ml of BJM RCS
- 10 Automix Syringe Mix Tips and Intraoral Tips
- 1 Mixing Pad

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1. Antibacterial mechanism of novel endodontic sealer, D. Kesler Shvero, N. Zaltsman, E. Weiss, N. Beyth, Hadassah School of Dental Medicine, Hebrew University, IADR Israeli Division Meeting, Tel-Aviv, June 2013.
2. Root canal sealers as Biofilm prevention: facts and speculations, M. Solomonov, **Эндодонтия**, Том VII, No. 1-2, 2014.
3. Antibiofilm Activity of Epoxy Sealer with Quaternary Ammonium Macromolecule, T. Becker, M. Solomonov, N. Sterer, R. Bar-Ness, A. Levin, A. Shemesh, The Maurice and Gabriela Goldschleger School of Dental Medicine Tel Aviv University, Program Number 0219, PER-IADR Congress, Jerusalem, Israel, 2016.
4. Dr. Michael Solomonov, Clinical Cases Report, March 2017.
5. Evaluating the physical properties of one novel and two well-established epoxy resin-based root canal sealers, M. Solomonov, J.B. Itzhak, Quintessence Publishing Endo 2017; 11 (4): 285-290.
6. Antibiofilm activity of epoxy sealer incorporated with quaternary ammonium macromolecule, T. Becker, N. Sterer, R. Bar-Ness, T. Toledano and M. Solomonov Evidence-Based Endodontics 2019; 4 (1): 1-6.
7. Dr. Gabriel Kaplan, Clinical Cases Report, March 2019.

For more state of the art products please
visit us at our new web site
www.bjmlabs.com

BJM LAB

B.J.M. Laboratories Ltd.
12 Hassadna St., Industrial Park,
Or-Yehuda 6022011, Israel
tel. 972-74-7000111 fax. 972-3-7353020
Info@bjmlabs.com www.bjmlabs.com

QUASER S.r.l.
Via Pian due Torri 19, 00146 Rome, Italy
tel. 39-06-55282325 fax. 39-06-55290821
info@quasersrl.it www.quasersrl.it