

BJM Q-Match *Smart Shade Matching Universal Restorative Composite*

ADVANTAGES

Smart
blending
properties

Simplified 3
shades
selection

Superb
handling

Universal

Excellent
strength
and
durability

FORMULA MAIN BLOCKS

Resin

Fillers

Photo-
initiators

Additives

Bis-GMA

Dental Glass

CQ

Pigments

TEGDMA

Fumed Silica

TPO

EDB

TESTS RESULTS*

**BJM
LAB**

3M

Filtek Suprem
Ultra



Tokuyama
Omnichroma



Dentsply
Sirona

TPH Spectra ST

Kerr

SimpliShade

Depth of Cure (mm)

9

9

9

9

9

Flexural strength (MPa)

106 ± 10

78 ± 10

73 ± 6

105 ± 10

79 ± 10

Compressive Strength (MPa)

289 ± 17

278 ± 28

315 ± 27

270 ± 27

260 ± 23

SBS** to Enamel (MPa)

27 ± 3

32 ± 3

24 ± 3

25 ± 2

22 ± 4

Water sorption (µg/mm³)

12 ± 1

20 ± 3

17 ± 2

15 ± 2

19 ± 2

Water solubility (µg/mm³)

1 ± 0.1

3 ± 2

1 ± 0.2

3 ± 1

2 ± 0.3

Radiopacity (mm Al)

2.74

2.14

1.27

2.43

2.50

Polymerization Shrinkage (%)

1.3

1.1

1.3

2.1

3.8

Degree of Conversion (%)

60

41

46

41

41

*All results and findings presented in this document are based on internal testing conducted by BJM LAB. These results may vary under different test conditions and should be interpreted accordingly.

** Using 3M Scotchbond universal adhesive

BJM LAB