


ADVANTAGES



FORMULA MAIN BLOCKS



HEMA

UDMA



Ethanol

Water



CQ

BAPO

EDB



Silane

Adhesion promoter

10 MDP

Filler

BONDING TEST RESULTS



Enamel eMAX CAD
Dentin 3M Lava
Zirconia



Composite
Core Buildup
Adhesive Resin Cement

	BJM LAB	Dentsply Sirona	3M
*w/o etching **Sandblasting			
Enamel*	18 ± 2	20 ± 3	16 ± 3
Dentin*	28 ± 3	25 ± 7	23 ± 8
Zirconia*	23 ± 2	18 ± 2	21 ± 4
eMAX CAD**	23 ± 2	15 ± 2	20 ± 2
3M Lava**	25 ± 3	18 ± 2	17 ± 2

*All results and findings presented in this document are based on internal testing conducted by BJM LAB. These results may vary under different test conditions and should be interpreted accordingly.